


## Add segmentation group

To add a segmentation group, follow these steps:



1. Click , the following window appears:

The screenshot shows a dialog box titled "New Segmentation Group" with a close button (X) in the top right corner. The dialog contains two input fields: "Segmentation Group\*" with a character count of "0 / 100" and "Segmentations". At the bottom, there are three buttons: "CANCEL", "SAVE", and "SAVE AND NEW".

2. Enter the segmentation group name,
3. Select the segmentations,
4. Click **Save** to add the group and close the window, or **Save and New** to add it and open a new window.