


## Assign segmentation to article

To assign a segmentation to an article, follow these steps:

1. Click , the following window appears:



The screenshot shows a 'New Record' dialog box. The title bar is blue with the text 'New Record' and a close button (X). The main area is white and contains a text input field with the placeholder text 'Segmentation Id \*'. At the bottom of the dialog, there are three buttons: 'CANCEL' (dark blue), 'SAVE' (light gray), and 'SAVE AND NEW' (light gray).

You can also assign a segmentation by selecting “Create Segmentation” from the Actions menu next to article.

2. Click the “Segmentation Id” field then select the segmentation you want,
3. Click **Save** to add the segmentation and close the window, or **Save and New** to add it and open a new window. The segmentation is added to the list.

You can also assign a segmentation to multiple articles as mentioned in “[Set Segmentation to Articles](#)”.

### Note:

- This button is not available if the selected catalog is a child catalog with the “Segmentation Mirroring” option activated.